

# Solido



Patient Overbed Table - Server

# Solido

The Solido, patient overbed table has a lightweight design and features an eating board with height adjustment and tilting.



- HEIGHT ADJUSTMENT
- TWO PIECE TILTING EATING BOARD
- SOLID AND STABLE CONSTRUCTION



## UNIVERSAL USE

- The ease of manipulation and ideal size of the eating board allows it to be used as a stable surface by the patient when eating.
- Two-piece tilting eating board option allows the patient to position the table at an angle for activities such as reading or writing.
- Flexibility of the Solido ensures it can be used by the patient both in bed and when seated in a chair.
- Low height of undercarriage allows use of the Solido together with a wide range of bed frames.

## STABILITY AND ROBUSTNESS

- Solid construction ensures comfortable and stable use for when eating or reading.
- Thanks to the solid frame it is possible to load the eating board with up to 30kg in the centre.
- Raised edges on the eating board provide a preventative measure against items slipping off the surface.
- Smooth and non-absorbent surfaces ensure that the Solido is easy to clean and decontaminate.
- Double castors and lightweight design make the Solido easy to manoeuvre and control.



## TECHNICAL DATA

Eating board dimensions	87×43 cm
Height of eating board	70 – 103 cm
Underpassing height (castors diameter 7,5 cm)	12,8 cm
Underpassing height (castors diameter 5 cm)	10,8 cm
Max. load in the centre of eating board	30 kg